

**Product Name: Phospho-Vimentin (Ser72) Rabbit
Monoclonal antibody
Catalog #: AMRe02762**

Summary

Production Name	Phospho-Vimentin (Ser72) Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	VIM
Alternative Names	VIM; Vimentin
Gene ID	7431
SwissProt ID	P08670.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Background


**Product Name: Phospho-Vimentin (Ser72) Rabbit
Monoclonal antibody
Catalog #: AMRe02762**

Vimentin an intermediate filament protein. Intermediate filament proteins are expressed in a tissue-specific manner. Desmin is the subunit specific for muscle and vimentin the subunit specific for mesenchymal tissue.

Research Area

Neuroscience

Image Data



Western blot analysis of Vimentin (Phospho-Ser72) in rat Brain lysates using Phospho-Vimentin (Ser72) antibody.

Note

For research use only.