

Product Name: Wnt2b Rabbit Monoclonal Antibody
Catalog #: AMRe02774



Summary

Production Name	Wnt2b Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	WNT2B
Alternative Names	WNT13; Wnt2b; XWNT2
Gene ID	7482
SwissProt ID	Q93097.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Background

Product Name: Wnt2b Rabbit Monoclonal Antibody
Catalog #: AMRe02774




Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.


Research Area

Stem Cells

Image Data



Immunocytochemistry analysis of Wnt2b (green) in 293t using Wnt2b antibody, and DAPI (blue).



Western blot analysis of Wnt2b in K562, C6, 3T3, HeLa lysates using Wnt2b antibody.

Note

For research use only.