

Summary

Production Name	IRF7 Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	IRF7
Alternative Names	IRF7; Interferon regulatory factor 7; IRF-7; IRF7A; IRF-7H
Gene ID	3665
SwissProt ID	Q92985.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Background

Product Name: IRF7 Rabbit Monoclonal Antibody
Catalog #: AMRe02793




Binds to the Q promoter (Qp) of EBV nuclear antigen 1 a (EBNA1) and may play a role in the regulation of EBV latency. Can activate distinct gene expression programs in macrophages and regulate the anti-tumor properties of primary macrophages.


Research Area

Immunology


Image Data



Immunocytochemistry analysis of IRF7 (green) in 293 using IRF7 antibody, and DAPI (blue).



Western blot analysis of IRF7 in rat Brain lysates using IRF7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using IRF7 antibody. High-pressure and

Product Name: IRF7 Rabbit Monoclonal Antibody
Catalog #: AMRe02793



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.