

Summary

Production Name	MST1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	MST1
Alternative Names	HGFL; MSP; MST1; NF15S2
Gene ID	4485
SwissProt ID	P26927.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 80 kDa; Observed MW: 96 kDa

Background

Product Name: MST1 Rabbit Monoclonal Antibody
Catalog #: AMRe02805




Probably has no proteolytic activity, since crucial characteristic of serine proteases catalytic sites are not conserved.


Research Area

Cell Biology


Image Data



Western blot analysis of MST1 in C2C12 lysates using MST1 antibody.



Western blot analysis of MST1 in HeLa, C6, rat Brain lysates using MST1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MST1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: MST1 Rabbit Monoclonal Antibody
Catalog #: AMRe02805



Note

For research use only.