

Product Name: ETS1 Rabbit Monoclonal Antibody**Catalog #: AMRe02807**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

Antigen Information

Gene Name	ETS1
Alternative Names	ETS1; EWSR2; Protein C-ets-1; p54
Gene ID	2113
SwissProt ID	P14921
Immunogen	A synthetic peptide of human ETS1


Background

Transcription factor.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ETS1 in HeLa, A549, HL-60, U2OS lysates using ETS1 antibody.