

Product Name: BDKRB2 Rabbit Monoclonal Antibody
Catalog #: AMRe02808



Summary

Production Name	BDKRB2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	BDKRB2
Alternative Names	B2BKR; B2BRA; B2R; BDKRB2; BK2 receptor; BK2R; BKR2; BR B2; BRB2; Kinin B2
Gene ID	624
SwissProt ID	P30411.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 80-100 kDa

Background

Receptor for bradykinin. It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger

Product Name: BDKRB2 Rabbit Monoclonal Antibody
Catalog #: AMRe02808




system.

Research Area

Cardiovascular

Image Data



Western blot analysis of BDKRB2 in A549, HL-60, U2OS lysates using BDKRB2 antibody.

Note

For research use only.