
Product Name: OR1A1 Rabbit Monoclonal Antibody**Catalog #: AMRe02814**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.65mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Antigen Information

Gene Name	OR1A1
Alternative Names	OR17-7
Gene ID	8383
SwissProt ID	Q9P1Q5
Immunogen	A synthetic peptide of human OR1A1


Background

Odorant receptor.


Research Area

Neuroscience

Image Data



Western blot analysis of OR1A1 in 3T3, MEF lysates using OR1A1 antibody.



Western blot analysis of OR1A1 in CHO-K1, C6, Jurkat lysates using OR1A1 antibody.