
Product Name: RBM43 Rabbit Monoclonal Antibody**Catalog #: AMRe02817**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Antigen Information


Gene Name	RBM43
Alternative Names	C2orf38
Gene ID	375287
SwissProt ID	Q6ZSC3
Immunogen	Recombinant protein of human RBM43

Background

Predicted to enable RNA binding activity.

Research Area

Image Data



Western blot analysis of RBM43 in HeLa lysates using RBM43 antibody.