

Summary

Production Name	FGFR3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FGFR3
Alternative Names	FGFR3; JTK4; Fibroblast growth factor receptor 3; FGFR-3; CD antigen CD333
Gene ID	2261
SwissProt ID	P22607.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 125 kDa

Background

Product Name: FGFR3 Rabbit Monoclonal Antibody
Catalog #: AMRe02820




Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.


Research Area

Cardiovascular


Image Data



Western blot analysis of FGFR3 in HeLa lysates using FGFR3 antibody.



Immunocytochemistry analysis of FGFR3 (green) in HeLa using FGFR3 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using FGFR3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: FGFR3 Rabbit Monoclonal Antibody
Catalog #: AMRe02820



Note

For research use only.