

Product Name: CAP2 Rabbit Monoclonal Antibody**Catalog #: AMRe02821**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Antigen Information

Gene Name	CAP2
Alternative Names	CAP 2; Adenylyl cyclase-associated protein 2
Gene ID	10486
SwissProt ID	P40123
Immunogen	A synthetic peptide of human CAP2


Background

May have a regulatory bifunctional role.

Research Area

Signal Transduction

Image Data



Western blot analysis of CAP2 in HeLa, CHO-K1 lysates using CAP2 antibody.