

Product Name: PACSIN3 Rabbit Monoclonal Antibody**Catalog #: AMRe02825**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 48 kDa; Observed MW: 48 kDa

Antigen Information


Gene Name	PACSIN3
Alternative Names	SDPIII
Gene ID	29763
SwissProt ID	Q9UKS6
Immunogen	A synthetic peptide of human PACSIN3

Background


Plays a role in endocytosis and regulates internalization of plasma membrane proteins.

Research Area


Image Data



Western blot analysis of PACSIN3 in HeLa, A549 lysates using PACSIN3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using PACSIN3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PACSIN3 (green) in HeLa using PACSIN3 antibody, and DAPI (blue)