
Product Name: CDK16 Rabbit Monoclonal Antibody**Catalog #: AMRe02832**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 65 kDa

Antigen Information

Gene Name	CDK16
Alternative Names	PCTK1; PCTAIRE; PCTAIRE1; PCTGAIRE
Gene ID	5127
SwissProt ID	Q00536
Immunogen	Recombinant protein of human PCTAIRE1


Background

Protein kinase that plays a role in vesicle-mediated transport processes and exocytosis.

Research Area

Signal Transduction

Image Data



Western blot analysis of PCTAIRE1 in A549, HL-60, U251 lysates using CDK16 antibody.