

**Product Name: Phospho-AKT1 (Thr450) Rabbit
Monoclonal Antibody
Catalog #: AMRe02836**

Summary

Production Name	Phospho-AKT1 (Thr450) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	AKT1
Alternative Names	AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha
Gene ID	207
SwissProt ID	P31749.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 56 kDa

**Product Name: Phospho-AKT1 (Thr450) Rabbit
Monoclonal Antibody
Catalog #: AMRe02836**


Background

Akt, also referred to as PKB or Rac, plays a critical role in controlling survival and apoptosis. This protein kinase is activated by insulin and various growth and survival factors to function in a wortmannin-sensitive pathway involving PI3 kinase. Akt is activated by phospholipid binding and activation loop phosphorylation at Thr308 by PDK1 and by phosphorylation within the carboxy terminus at Ser473.

Research Area

Neuroscience

Image Data



Western blot analysis of Phospho-AKT1 (Thr450) in U2OS, C6 lysates using Phospho-AKT1 (Thr450) antibody.

Note

For research use only.