

**Product Name: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**
Catalog #: AMRe02850

Summary

Production Name	Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	EZR
Alternative Names	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
Gene ID	7430
SwissProt ID	P15311.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 69 kDa; Observed MW: 75,80 kDa

**Product Name: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody
Catalog #: AMRe02850**




Background

Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.


Research Area

Cell Biology

Image Data




Western blot analysis of Phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) in 3T3, 3T3+PDGF lysates using Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) antibody.



Western blot analysis of Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) in HeLa, CHO-K1, C6 lysates using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody.

**Product Name: Phospho-Ezrin/Radixin/Moesin
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody
Catalog #: AMRe02850**



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.