

Product Name: Phospho-Filamin A/B (Ser2152/Ser2107)
Rabbit Monoclonal Antibody
Catalog #: AMRe02852

Summary

Production Name	Phospho-Filamin A/B (Ser2152/Ser2107) Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FLNB
Alternative Names	FLN; FMD; MNS; OPD; ABPX; CSBS; CVD1; FGS2; FLN1; NHBP; OPD1; OPD2; XLVD; XMVD; FLN-A; ABP-280
Gene ID	2317
SwissProt ID	O75369.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 278 kDa; Observed MW: 281 kDa


Product Name: Phospho-Filamin A/B (Ser2152/Ser2107)
Rabbit Monoclonal Antibody
Catalog #: AMRe02852

Background


Connects cell membrane constituents to the actin cytoskeleton.

Research Area

Image Data




Immunocytochemistry analysis of Phospho-Filamin A/B (Ser2152/Ser2107) (green) in A549 using Phospho-Filamin A/B (Ser2152/Ser2107) antibody, and DAPI (blue).



Western blot analysis of Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) in HeLa, A549, U2OS, C6 lysates using Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) antibody.

**Product Name: Phospho-Filamin A/B (Ser2152/Ser2107)
Rabbit Monoclonal Antibody
Catalog #: AMRe02852**



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.