

Product Name: Bora Rabbit Monoclonal Antibody**Catalog #: AMRe02917**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 61 kDa; Observed MW: 61 kDa

Antigen Information

Gene Name	BORA
Alternative Names	C13orf34
Gene ID	79866
SwissProt ID	Q6PGQ7
Immunogen	Recombinant protein of human Bora


Background

Required for the activation of AURKA at the onset of mitosis.

Research Area

Cell Biology

Image Data



Western blot analysis of Bora in HL-60, C6 lysates using Bora antibody.