

Product Name: Cadherin 10 Rabbit Monoclonal Antibody
Catalog #: AMRe02921

Summary

Production Name	Cadherin 10 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	CDH10
Alternative Names	T2-cadherin; CDH10
Gene ID	1008
SwissProt ID	Q9Y6N8.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 110 kDa

Product Name: Cadherin 10 Rabbit Monoclonal Antibody
Catalog #: AMRe02921




Background

Cadherins are calcium-dependent cell adhesion proteins.


Research Area

Cell Biology

Image Data



Western blot analysis of cadherin 10 in MEF lysates using Cadherin 10 antibody.



Western blot analysis of cadherin 10 in U2OS, C6 lysates using cadherin 10 antibody.

Note

For research use only.