

Product Name: Integrin alpha L Rabbit Monoclonal Antibody**Catalog #: AMRe02932**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.64mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 128 kDa; Observed MW: 160 kDa

Antigen Information

Gene Name	Itgal
Alternative Names	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Gene ID	16408.0
SwissProt ID	P24063
Immunogen	Recombinant protein of mouse CD11a


Background

Integrin ITGAL/ITGB2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4.


Research Area

Tags & Cell Markers

Image Data



Immunohistochemistry analysis of paraffin-embedded mouse lung using Integrin alpha L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD11a in A549, HL-60 lysates using CD11a antibody.