
Product Name: CD229 Rabbit Monoclonal antibody**Catalog #: AMRe02933**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.18mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 72 kDa; Observed MW: 72 kDa

Antigen Information


Gene Name	LY9
Alternative Names	hly9; mLY9; CD229; SLAMF3
Gene ID	4063
SwissProt ID	Q9HBG7
Immunogen	Recombinant protein of human CD229

Background

Self-ligand receptor of the signaling lymphocytic activation molecule (SLAM) family.

Research Area

Image Data



Western blot analysis of CD229 in HeLa, A549, HL-60, U251 lysates using CD229 antibody.