

Product Name: Cellubrevin Rabbit Monoclonal Antibody**Catalog #: AMRe02941**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.63mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

Antigen Information

Gene Name	VAMP3
Alternative Names	CEB; Cellubrevin; Synaptobrevin 3; VAMP3
Gene ID	9341
SwissProt ID	Q15836
Immunogen	Recombinant protein of human Cellubrevin


Background

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.


Research Area

Signal Transduction

Image Data



Western blot analysis of ubrevin in HeLa, A549, MCF-7, 293T lysates using ubrevin antibody.



Immunocytochemistry analysis of Cellubrevin (green) in hela using Cellubrevin antibody, and DAPI (blue)