

Product Name: Clathrin Heavy Chain Rabbit Monoclonal antibody
Catalog #: AMRe02945

Summary

Production Name	Clathrin Heavy Chain Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	CLTC
Alternative Names	Hc; CHC; CHC17; MRD56; CLH-17; CLTCL2
Gene ID	1213
SwissProt ID	Q00610.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 192 kDa; Observed MW: 192 kDa

Product Name: Clathrin Heavy Chain Rabbit Monoclonal antibody
Catalog #: AMRe02945


Background

Clathrin is the major protein of the polyhedral coat of coated pits and vesicles.


Research Area

Signal Transduction

Image Data



Western blot analysis of Clathrin heavy chain in Hela, A549, HL-60, U251, U87-MG lysates using Clathrin Heavy Chain antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Clathrin heavy chain antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.