
Product Name: CSDE1 Rabbit Monoclonal antibody**Catalog #: AMRe02954**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.6mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 89 kDa; Observed MW: 89 kDa

Antigen Information

Gene Name	CSDE1
Alternative Names	UNR; D1S155E
Gene ID	7812
SwissProt ID	O75534
Immunogen	A synthetic peptide of human CSDE1/NRU


Background

RNA-binding protein involved in translationally coupled mRNA turnover.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of CSDE1/NRU in A549, HL-60 lysates using CSDE1 antibody.