

Product Name: ERCC8 Rabbit Monoclonal antibody
Catalog #: AMRe02990



Summary

Production Name	ERCC8 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	ERCC8
Alternative Names	CKN1; CSA; ERCC8
Gene ID	1161
SwissProt ID	Q13216.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Background

Substrate-recognition component of the CSA complex, a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex,

Product Name: ERCC8 Rabbit Monoclonal antibody
Catalog #: AMRe02990




involved in transcription-coupled nucleotide excision repair.

Research Area

-

Image Data



Western blot analysis of ERCC8 in A549, HL-60, U251, U87-MG lysates using ERCC8 antibody.

Note

For research use only.