
Product Name: EYA2 Rabbit Monoclonal Antibody**Catalog #: AMRe02991**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 59 kDa; Observed MW: 59 kDa

Antigen Information

Gene Name	EYA2
Alternative Names	EAB1
Gene ID	2139
SwissProt ID	O00167
Immunogen	A synthetic peptide of human EYA2


Background

Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.