

Product Name: FHL1 Rabbit Monoclonal Antibody
Catalog #: AMRe02999



Summary

Production Name	FHL1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FHL1
Alternative Names	FHL1; FHL1A; FHL1B; FLH1A; RBMX1A; RBMX1B; SLIM1; SLIMMER; XMPMA
Gene ID	2273
SwissProt ID	Q13642.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 36 kDa

Background

Product Name: FHL1 Rabbit Monoclonal Antibody
Catalog #: AMRe02999




May have an involvement in muscle development or hypertrophy.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of FHL1 in HeLa, U251 lysates using FHL1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FHL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.