

Summary

Production Name	FKBP52 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FKBP4
Alternative Names	HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PPlase
Gene ID	2288
SwissProt ID	Q02790.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

Product Name: FKBP52 Rabbit Monoclonal Antibody
Catalog #: AMRe03000




Immunophilin protein with PPlase and co-chaperone activities (By similarity). Component of unliganded steroid receptors heterocomplexes through interaction with heat-shock protein 90 (HSP90).


Research Area

Signal Transduction


Image Data



Western blot analysis of FKBP52 in Raw264.7 lysates using FKBP52 antibody.



Western blot analysis of FKBP52 in U2OS, C6 lysates using FKBP52 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FKBP52 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: FKBP52 Rabbit Monoclonal Antibody
Catalog #: AMRe03000



Note

For research use only.