

Product Name: Flotillin 1 Rabbit Monoclonal Antibody**Catalog #: AMRe03002**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Antigen Information

Gene Name	FLOT1
Alternative Names	Flotillin-1
Gene ID	10211
SwissProt ID	O75955
Immunogen	A synthetic peptide of human Flotillin 1


Background

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Research Area

Cell Biology

Image Data



Western blot analysis of Flotillin 1 in HeLa, A549, MCF-7, U251 lysates using Flotillin 1 antibody.