

Product Name: GMEB1 Rabbit Monoclonal Antibody
Catalog #: AMRe03013



Summary

Production Name	GMEB1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GMEB1
Alternative Names	PIF96; P96PIF
Gene ID	10691
SwissProt ID	Q9Y692.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 63 kDa; Observed MW: 80 kDa

Background

Trans-acting factor that binds to glucocorticoid modulatory elements (GME) present in the TAT (tyrosine aminotransferase)

Product Name: GMEB1 Rabbit Monoclonal Antibody
Catalog #: AMRe03013




promoter and increases sensitivity to low concentrations of glucocorticoids.

Research Area

-

Image Data



Western blot analysis of GMEB1 in A549, HL-60, U2OS lysates using GMEB1 antibody.

Note

For research use only.