

Summary

Production Name	GRSF1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GRSF1
Alternative Names	G-rich sequence factor 1
Gene ID	2926
SwissProt ID	Q12849.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Background

Product Name: GRSF1 Rabbit Monoclonal Antibody
Catalog #: AMRe03016




Regulator of post-transcriptional mitochondrial gene expression, required for assembly of the mitochondrial ribosome and for recruitment of mRNA and lncRNA.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Immunocytochemistry analysis of GRSF1 (green) in CEM using GRSF1 antibody, and DAPI (blue).



Western blot analysis of GRSF1 in HeLa, Ramos lysates using GRSF1 antibody.

Note

For research use only.