

Product Name: GTPBP4 Rabbit Monoclonal antibody
Catalog #: AMRe03017



Summary

Production Name	GTPBP4 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	GTPBP4
Alternative Names	CRFG; GTP binding protein 4; Gtpbp4; NGB; NOG1
Gene ID	23560
SwissProt ID	Q9BZE4.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 74 kDa; Observed MW: 80 kDa

Background

Product Name: GTPBP4 Rabbit Monoclonal antibody
Catalog #: AMRe03017




Involved in the biogenesis of the 60S ribosomal subunit.


Research Area

Signal Transduction

Image Data



Western blot analysis of GTPBP4/NOG1 in A549, HL-60 lysates using GTPBP4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using GTPBP4/NOG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.