

Product Name: HDLBP Rabbit Monoclonal Antibody
Catalog #: AMRe03021



Summary

Production Name	HDLBP Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	HDLBP
Alternative Names	HBP; VGL; PRO2900
Gene ID	3069
SwissProt ID	Q00341.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 141 kDa; Observed MW: 141 kDa

Background

Appears to play a role in cell sterol metabolism.


Product Name: HDLBP Rabbit Monoclonal Antibody
Catalog #: AMRe03021



Research Area

Cardiovascular

Image Data



Western blot analysis of HDLBP in HeLa, A549, HL-60, U251, U87-MG lysates using HDLBP antibody.

Note

For research use only.