

Product Name: HJURP Rabbit Monoclonal Antibody
Catalog #: AMRe03030



Summary

Production Name	HJURP Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	HJURP
Alternative Names	FAKTS; URLC9; hFLEG1
Gene ID	55355
SwissProt ID	Q8NCD3.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 84 kDa

Background


Product Name: HJURP Rabbit Monoclonal Antibody
Catalog #: AMRe03030




Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres.

Research Area

Image Data



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.