
Product Name: LILRB1 Rabbit Monoclonal Antibody**Catalog #: AMRe03066**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 71 kDa; Observed MW: 110 kDa

Antigen Information

Gene Name	LILRB1
Alternative Names	ILT2; LIR1; MIR7; PIRB; CD85J; ILT-2; LIR-1; MIR-7; PIR-B
Gene ID	10859
SwissProt ID	Q8NHL6
Immunogen	A synthetic peptide of human LILRB1


Background

Receptor for class I MHC antigens.

Research Area

Immunology

Image Data



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.