
Product Name: LINE 1 ORF1p Rabbit Monoclonal Antibody**Catalog #: AMRe03068**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Antigen Information

Gene Name	L1RE1
Alternative Names	L1ORF1p; L1RE1; LRE1
Gene ID	105370531.0
SwissProt ID	Q9UN81
Immunogen	


Background

Nucleic acid-binding protein which is essential for retrotransposition of LINE-1 elements in the genome.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of LINE1 ORF1p in HeLa, A549 lysates using LINE 1 ORF1p antibody.