

Product Name: LLGL1 Rabbit Monoclonal Antibody
Catalog #: AMRe03069



Summary

Production Name	LLGL1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	LLGL1
Alternative Names	DLG4; HUGL; LLGL; Lgl1; Mgl1; HUGL1; HUGL-1
Gene ID	3996
SwissProt ID	Q15334.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 115 kDa; Observed MW: 115 kDa

Background

Cortical cytoskeleton protein found in a complex involved in maintaining cell polarity and epithelial integrity.


Product Name: LLGL1 Rabbit Monoclonal Antibody
Catalog #: AMRe03069



Research Area

Neuroscience

Image Data



Western blot analysis of LLGL1 in A549, U2OS lysates using LLGL1 antibody.

Note

For research use only.