

Product Name: MAPKAP Kinase 3 Rabbit Monoclonal antibody
Catalog #: AMRe03088



Summary

Production Name	MAPKAP Kinase 3 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	MAPKAPK3
Alternative Names	MAPKAPK3; MAP kinase-activated protein kinase 3; MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAP-K3; MAPKAPK-3; MK-3; Chromosome 3p kinase; 3pK
Gene ID	7867
SwissProt ID	Q16644.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Product Name: MAPKAP Kinase 3 Rabbit Monoclonal antibody
Catalog #: AMRe03088




Background

Stress-activated serine/threonine-protein kinase involved in cytokines production, endocytosis, cell migration, chromatin remodeling and transcriptional regulation.

Research Area

Signal Transduction

Image Data



Western blot analysis of MK3 in A549, HL-60 lysates using MAPKAP Kinase 3 antibody.

Note

For research use only.