

Product Name: MLLT6 Rabbit Monoclonal Antibody
Catalog #: AMRe03090



Summary

Production Name	MLLT6 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	MLLT6
Alternative Names	AF17
Gene ID	4302
SwissProt ID	P55198.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 112 kDa; Observed MW: 112 kDa

Background

Enables histone binding activity and nucleosome binding activity. Predicted to be involved in positive regulation of

Product Name: MLLT6 Rabbit Monoclonal Antibody
Catalog #: AMRe03090




transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including negative regulation of histone H3-K79 methylation; renal potassium excretion; and renal sodium excretion. Predicted to be active in nucleus.

Research Area

-

Image Data



Western blot analysis of MLLT6 in HeLa, A549, HL-60, U2OS lysates using MLLT6 antibody.

Note

For research use only.