

Product Name: MPP1 Rabbit Monoclonal Antibody**Catalog #: AMRe03092**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 55 kDa

Antigen Information


Gene Name	MPP1
Alternative Names	MRG1; PEMP; AAG12; EMP55; DXS552E
Gene ID	4354
SwissProt ID	Q00013
Immunogen	A synthetic peptide of human MPP1

Background

Essential regulator of neutrophil polarity.

Research Area

Image Data



Western blot analysis of MPP1 in A549, HL-60 lysates using MPP1 antibody.