

**Product Name: Interferon Regulatory Factor 4 Rabbit
Monoclonal Antibody**
Catalog #: AMRe03098



Summary

Production Name	Interferon Regulatory Factor 4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	IRF4
Alternative Names	IRF4; MUM1; Interferon regulatory factor 4; IRF-4; Lymphocyte-specific interferon regulatory factor; LSIRF; Multiple myeloma oncogene 1; NF-EM5
Gene ID	3662
SwissProt ID	Q15306.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

**Product Name: Interferon Regulatory Factor 4 Rabbit
Monoclonal Antibody
Catalog #: AMRe03098**




Background

Melanoma associated antigen (mutated) 1 (MUM1, EXPAND1) is a PWWP-domain containing chromatin binding protein involved in maintaining chromatin architecture of interphase chromosomes.

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of MUM1 in Ramos lysates using Interferon Regulatory Factor 4 antibody.

Note

For research use only.