

Product Name: Nectin 2 Rabbit Monoclonal Antibody
Catalog #: AMRe03106



Summary

Production Name	Nectin 2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	NECTIN2
Alternative Names	HVEB; PRR2; CD112; PVRR2
Gene ID	5819
SwissProt ID	Q92692.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 70-80 kDa

Background

Probable cell adhesion protein.


Product Name: Nectin 2 Rabbit Monoclonal Antibody
Catalog #: AMRe03106



Research Area

Immunology

Image Data



Western blot analysis of Nectin 2 in Hela, A549, U2OS, U87-MG lysates using Nectin 2 antibody.

Note

For research use only.