

Product Name: NUP107 Rabbit Monoclonal Antibody
Catalog #: AMRe03122



Summary

Production Name	NUP107 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	NUP107
Alternative Names	ODG6; NUP84; NPHS11; ODG6; GAMOS7
Gene ID	57122
SwissProt ID	P57740.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 106 kDa; Observed MW: 106 kDa

Background

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance (PubMed:12552102, PubMed:15229283).

Product Name: NUP107 Rabbit Monoclonal Antibody
Catalog #: AMRe03122




Required for the assembly of peripheral proteins into the NPC (PubMed:15229283, PubMed:12552102). May anchor NUP62 to the NPC (PubMed:15229283).

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Nup107 in A549, HL-60, U251 lysates using NUP107 antibody.

Note

For research use only.