

Product Name: NUP50 Rabbit Monoclonal antibody
Catalog #: AMRe03124



Summary

Production Name	NUP50 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	NUP50
Alternative Names	NPAP60; NPAP60L
Gene ID	10762
SwissProt ID	Q9UKX7.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

Background

Product Name: NUP50 Rabbit Monoclonal antibody
Catalog #: AMRe03124




Component of the nuclear pore complex that has a direct role in nuclear protein import.


Research Area

-

Image Data



Immunocytochemistry analysis of NUP50 (green) in Jurkat using NUP50 antibody, and DAPI (blue).



Western blot analysis of NUP50 in HeLa, A549, HL-60, U251 lysates using NUP50 antibody.

Note

For research use only.