

Product Name: NUP98 Rabbit Monoclonal antibody
Catalog #: AMRe03126



Summary

Production Name	NUP98 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	NUP98
Alternative Names	96 kDa nucleoporin; 98 kDa nucleoporin; ADAR2; ADIR2; GLFG-repeat containing nucleoporin; Nuclear pore complex protein Nup96; Nuclear pore complex protein Nup98 Nup96; Nucleoporin 98kD; nucleoporin 98kDa; Nucleoporin Nup96; Nucleoporin Nup98; NUP196; NUP96; Nup98; Nup98-Nup96; NUP98_HUMAN.
Gene ID	4928
SwissProt ID	P52948.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 198 kDa; Observed MW: 98 kDa

Product Name: NUP98 Rabbit Monoclonal antibody
Catalog #: AMRe03126




Background

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are involved in the bidirectional transport across the NPC. May anchor NUP153 and TPR to the NPC.

Research Area

Signal Transduction

Image Data



Western blot analysis of NUP98 in HeLa, A549, HL-60 lysates using NUP98 antibody.

Note

For research use only.