

Summary

Production Name	PF4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PF4
Alternative Names	C-X-C motif chemokine 4; CXCL4; Iroplact; OncostatinA; PF4; Platelet factor 4; SCYB4; short form; Small inducible cytokine subfamily member 4
Gene ID	5196
SwissProt ID	P02776.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

Product Name: PF4 Rabbit Monoclonal Antibody
Catalog #: AMRe03147




Background

Released during platelet aggregation. Neutralizes the anticoagulant effect of heparin because it binds more strongly to heparin than to the chondroitin-4-sulfate chains of the carrier molecule. Chemotactic for neutrophils and monocytes. Inhibits endothelial cell proliferation, the short form is a more potent inhibitor than the longer form.


Research Area

Immunology


Image Data



Western blot analysis of PF4 in rat spleen lysates using PF4 antibody.



Western blot analysis of PF4 in A549, HL-60 lysates using PF4 antibody.



Product Name: PF4 Rabbit Monoclonal Antibody
Catalog #: AMRe03147



Immunohistochemistry analysis of paraffin-embedded liver cancer using PF4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.