

Product Name: PKC theta Rabbit Monoclonal Antibody
Catalog #: AMRe03151



Summary

Production Name	PKC theta Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PRKCQ
Alternative Names	PRKCQ; PRKCT; Protein kinase C theta type; nPKC-theta
Gene ID	5588
SwissProt ID	Q04759.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 82 kDa; Observed MW: 79 kDa

Background

This is a calcium-independent, phospholipid-dependent, serine- and threonine-specific enzyme. Essential for T-cell receptor

Product Name: PKC theta Rabbit Monoclonal Antibody
Catalog #: AMRe03151




(TCR)-mediated T-cell activation, but is dispensable during TCR-dependent thymocyte development.

Research Area

Signal Transduction

Image Data



Western blot analysis of PKC theta/PRKCQ in U2OS, C6 lysates using PKC theta antibody.

Note

For research use only.