
Product Name: Pleckstrin Rabbit Monoclonal antibody**Catalog #: AMRe03152**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Antigen Information

Gene Name	PLEK
Alternative Names	P47; Platelet 47 kDa protein
Gene ID	5341
SwissProt ID	P08567
Immunogen	A synthetic peptide of human Pleckstrin


Background

Major protein kinase C substrate of platelets.

Research Area

Cardiovascular

Image Data



Western blot analysis of Pleckstrin in HeLa, HL-60 lysates using Pleckstrin antibody.