

Product Name: PMF 1 Rabbit Monoclonal antibody
Catalog #: AMRe03153



Summary

Production Name	PMF 1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	PMF1
Alternative Names	Polyamine-modulated factor 1
Gene ID	100527963
SwissProt ID	Q6P1K2.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Product Name: PMF 1 Rabbit Monoclonal antibody
Catalog #: AMRe03153




Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis.


Research Area

-

Image Data



Western blot analysis of PMF1 in C6, 3T3, Raw264.7 lysates using PMF 1 antibody.



Western blot analysis of PMF1 in HeLa, A549 lysates using PMF1 antibody.

Note

For research use only.