

Product Name: Rab29 Rabbit Monoclonal Antibody
Catalog #: AMRe03174



Summary

Production Name	Rab29 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	RAB29
Alternative Names	RAB7L; RAB7L1
Gene ID	8934
SwissProt ID	O14966.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

The small GTPases Rab are key regulators in vesicle trafficking.


Product Name: Rab29 Rabbit Monoclonal Antibody
Catalog #: AMRe03174



Research Area

-

Image Data



Western blot analysis of RAB29 in HeLa, A549, HL-60, U251 lysates using Rab29 antibody.

Note

For research use only.