

Product Name: RGS6 Rabbit Monoclonal Antibody**Catalog #: AMRe03185**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Antigen Information


Gene Name	RGS6
Alternative Names	regulator of G protein signaling 6; GAP; S914; HA117
Gene ID	9628
SwissProt ID	P49758
Immunogen	A synthetic peptide of human RGS6

Background

Regulates G protein-coupled receptor signaling cascades.

Research Area

Image Data



Western blot analysis of RGS6 in A549, HL-60, U2OS lysates using RGS6 antibody.