

Product Name: SAP97 Rabbit Monoclonal antibody
Catalog #: AMRe03194



Summary

Production Name	SAP97 Rabbit Monoclonal antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	DLG1
Alternative Names	DLG1; DLGH1; hDlg; SAP97
Gene ID	1739
SwissProt ID	Q12959.

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 100 kDa; Observed MW: 100 kDa

Background

Product Name: SAP97 Rabbit Monoclonal antibody
Catalog #: AMRe03194




Essential multidomain scaffolding protein required for normal development (By similarity). Recruits channels, receptors and signaling molecules to discrete plasma membrane domains in polarized cells.


Research Area

Neuroscience

Image Data



Western blot analysis of SAP97 in mouse brain, rat brain lysates using SAP97 antibody.



Western blot analysis of SAP97 in HeLa, A549, HL-60 lysates using SAP97 antibody.

Note

For research use only.